

Global Cell Solutions

Adhesion Assay Kit

Introduction:

The Global Eukaryotic Microcarrier (GEM) is a new generation culture surface for transformed, primary, and stem cell culture. Composed of alginate and coated in adhesion molecules, the GEM transforms adherent cells into a pipette-able, magnetic substrate. Sample, dispense, and transfer some or all of your cells just like a reagent. With the CubeMagnet or BioMagnet you can hold, pull, levitate, agitate and even rotate cells like never before. Now combine that with the ability to transfect, cryopreserve, assay and capture fluorescent images without removing the cells and you truly have the next generation of cell culture.


This Adhesion Assay Kit includes all the materials needed to get started: GEM, a low adhesion plate, Accutase, and magnets. The Adhesion Assay Kit contains all the available GEM coatings for cell culture: gelatin, collagen I, collagen IV, laminin, basement membrane, fibronectin, and poly-D-lysine. Use these coated GEMs in place of traditional culture matrices.

Questions? Suggestions? We want to help! Please contact us at:
+1-434-975-4271, ext 103
support@globalcellsolutions.com

The 24-well ultra low adhesion (ULA) plate is specially prepared to inhibit the attachment of cells ensuring adhesion to the GEM. Additional sizes and layouts are available on our web site depending on your application.

Two **strong** magnets are included as well. The stronger Cube Magnet is ideal for manipulating GEMs in solution. Rotate, hold or pull down GEMs easily. For levitating GEM in the ULA plate, use the BioMagnet. Its flat shape is ideal for applications requiring a broader application of the magnetic field.

Lastly, we've included the dissociation agent Accutase. Although we recommend keeping cells on the GEM throughout culture, assay, and application, we understand some applications require a cell suspension. Accutase is a trypsin alternative that is gentler on cells and doesn't require quenching with serum. Need a single cell suspension for flow cytometry or similar application? We also offer Accumax, which has all the dissociation power of Accutase while inhibiting the formation of cell-cell clumps in solution. See "Related Products" below for ordering information.



Catalog Number: GKT-5001

Kit Contents:

- 0.5mL gelatin GEM
- 0.5mL collagen I GEM
- 0.5mL collagen IV GEM
- 0.5mL laminin GEM
- 0.5mL basement membrane GEM
- 0.5mL fibronectin GEM
- 0.5mL poly-D-lysine GEM
- One (1) 24-well Ultra Low Adhesion (ULA) plate
- One (1) CubeMagnet
- One (1) BioMagnet
- One (1) bottle Accutase (100 mL)

Note: Do NOT thaw Accutase at 37° C (e.g. in a water bath). Thaw at room temperature or overnight at 4° C.

Storage and Disposal:

- **GEM** Store at 4° C and dispose of with biohazard waste.
- **CubeMagnet** Store away from electronics and other devices affected by magnetic fields. A used plastic bottle makes a great container.
- **BioMagnet** Store away from electronics and other devices affected by magnetic fields.
- **Accutase** Store at -20° C until opened. Store at 4° after opening. Dispose of with biohazard waste.

Related Protocols:

- Culturing Cells Using the GEM
- Immunohistochemical Staining of Cells on the GEM

Related Products:

Item	Catalog Number	Includes
GEM Supply Kit - Gelatin	GKT-7000G	2 mL GEM, 24-well ULA plate
GEM Supply Kit - Laminin	GKT-7000L	2 mL GEM, 24-well ULA plate
GEM Supply Kit - poly-D-lysine	GKT-7000PDL	2 mL GEM, 24-well ULA plate
Accutase	GSP-0101-100	100 mL bottle
Accumax	GSP-0102-100	100 mL bottle
60 mm ULA petri plate	GAA-2000	5 plates
100 mm ULA petri plate	GAA-2001	5 plates
6-well ULA plate	GAA-2002	1 plate
24-well plate	GAA-2003	1 plate
96-well plate	GAA-2004	1 plate